


Cylinder pressure regulator FM53

spectron



FM53-L-300-1,5-DIN10-M-M-0-A-B-N2



FM53-R-200-10-CGA540-M-M-CE6-A-B-Ar


Product features

- 2-stage diaphragm cylinder pressure regulator
- Safe and easy gas withdrawal from single cylinders with non-corrosive gases and gas mixtures up to quality 6.0
- 6 ports for flexible and individual configuration
- Very high control accuracy
- Integrated relief valve
- Metal-to-metal seal to atmosphere
- Central filter
- Approved for use with oxygen
- New laboratory-style design
- Including O-ring cylinder connection for mounting without tools
- Suitable for ECD-applications
- Pressure regulator can be evacuated

Technical data

Inlet pressure P_1	max. 300 bar
Max. outlet pressure P_2	1,5 / 4 / 10 / 20 bar
Materials	
Body:	chrome plated brass
Diaphragm:	Hastelloy C276
Valve seat:	PA
Soft goods:	FPM
Filter:	Sinterbronze
Inlet	Cylinder connection
Other In- / outlets	1/4"-NPT female
Temperature range	-20°C to +60°C
Leakrate (toatmosphere):	$\leq 10^{-8}$ mbar l/s He
(via seat):	$\leq 10^{-6}$ mbar l/s He
Weight	1,7 kg



Flow curves FM53



Cylinder pressure regulator FM53

spectrolab

Dimensions FM53





Ordering information:

Cylinder pressure regulator FM53

In- / outlets

CM 1/4	1/8" / 1/4" / 1/2"-compression fitting (brass)
CM 6	3 / 6 / 8 / 10 / 12 mm - compression fitting (brass)
CE 1/4	1/8" / 1/4" / 1/2"-compression fitting (stainless steel)
CE 6	3 / 6 / 8 / 10 / 12 mm - compression fitting (stainless steel)
B	plug
A	relief valve - 1/4"-NPT-F
xxxx	detailed specification of cylinder connection (e.g. BS no.3)
K	contact pressure gauge
M	pressure gauge
0	1/4"-NPT-F (port only)



FM53 - L - 300 - 20 - 1 - 2 - 3 - 4 - 5 - 6 - H2

Inlet

L	left
R	right

Inlet pressure P_1

100	max. 100 bar
200	max. 200 bar
300	max. 300 bar

Outlet pressure P_2

1,5	max. 1,5 bar
4	max. 4 bar
10	max. 10 bar
20	max. 20 bar

Gas type

Please specify gas type!

For inlet pressure up to 50 bar and outlet pressure of max. 1,5 or 4 bar please use M51.